Produktgruppe Messbecher

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Hersteller

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Article description

Measuring cup (all models)  
Designation: "M"

Intended use

Product for administering liquid or powdered medicines.

Target/patient group

Patients themselves (adults), a third party, or a healthcare professional can administer drugs. Only adults who are able to read and understand the meaning of the graduation scale should handle the measuring cups. Patients with physical or cognitive limitations require third-party assistance to use the measuring cup to ensure correct use and dosage. If in doubt or in case of dosage problems, consult your doctor, pharmacist, or specialist.

Composition of the medical device

Polypropylene\*

\*Please refer to the product specifications for exact device composition details.

Dosing deviation

These devices are compliant with Chapter 2.9.27 of the European Pharmacopoeia.

Total dosing deviation of the measured volume greater than 2,00 ml is  $\pm 20\%$ .

Total dosing deviation of the measured volume greater than 1,00 ml and less than 2,00 ml is  $\pm 30\%$ .

Total dosing deviation of the measured volume from 0,50 ml to 1,00 ml is  $\pm 35\%$ .

The dosing deviation was tested using a liquid medication: alcoholic solution, color: brown, viscosity: 3.1 mPas [23°C]. In addition, separate validations must be carried out using the planned medication (liquid / powder) to check suitability: e.g., dosing/withdrawal accuracy, adhesion, readability, etc.

Instructions for use

1. Check that the cup is in good condition and the scale is readable before each use. Ensure sufficient lighting conditions to be able to read the scale.
2. Pour the volume into the measuring cup according to prescribed dosage, placing the cup on a flat surface without vibration and aligning the level of the drug with the scale. The specification for the dosage of the liquid or powder is made by the distributor of the total packaging material and not by the manufacturer.
3. Administer the entire dosed volume.

CAUTION: gradually tilt the cup during administration in order to adjust the flow allowing the patient to swallow properly.

4. Clean the cup with water as specified in the cleaning instructions.
5. Store the cup in its original box.
6. Store the measuring cup including medication out of children's reach.
7. The measuring cup may only be combined with my matching/approved cap.
8. The measuring cup must not be used to open/close the closure system.

Cleaning

Clean the cup before each use. The measuring cup should be cleaned with water without any additions. It must be fully dried before reuse.

Produktgruppe Messbecher

Warning and safety instructions

- Do not use the cup if you notice any deterioration in its constituent material or shape (cracked or missing material). There is a risk of skin injury.
- Do not allow children to use the measuring cup unsupervised to avoid incorrect dosing.
- Do not bite the measuring cup to avoid possible skin injuries.
- The measuring cup must be clean before use to prevent infection.
- The measuring cup should only be used for its intended purpose. Dosing accuracy affects the effectiveness of medicines. This can prevent them from working or cause a delay in their effectiveness, which can have serious consequences for your health.
- Do not use the measuring cup for multiple patients to prevent infection. 
- If you lose the cup, do not use any other type of cup. Dosing accuracy affects the effectiveness of medicines. This can prevent them from working or cause a delay in their effectiveness, which can have serious consequences for your health.
- Do not use these devices if you know you are allergic to polypropylene\*. Allergic reactions may occur.  
\*Please refer to the product specifications for exact device composition details.

Side effects

Heinlein Plastik-Technik GmbH is currently unaware of any adverse effects connected to the use of these devices. If any incidents occur during the use of this medical device, please report them immediately:

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If serious incidents occur (serious injuries, deaths), please report them immediately to us and the responsible state authority.

Storage

Store in the original containers, protected against humidity and UV-radiation, storage temperature +5 to +35°C, at least 48 hours before processing at ambient temperature (20-25°C). For that purpose the cartons have to be opened.

Lifespan

The shelf life is 5 years from the production date. The expiry date can be found on the bag label

Disposal of the medical device

In the absence of any contrary country-specific disposal regulations, this medical device may be disposed of with household waste.

Communication to the patient

- Information on use, dosage deviation, and warnings/safety instructions must be communicated to patients in an appropriate and adequate manner.
- Do not use the measuring cup for multiple patients.
- Measuring cups can be used until the end of the medication's shelf life, after which they must be disposed of.

Symbols

Manufacturer logo  
Heinlein Plastik-Technik GmbH